



SIGNA™ PET/MR AIR™*

The clarity you seek.
The comfort they deserve.

Introducing the all-new SIGNA™ PET/MR AIR™

Whether you're reading clinical images, researching diagnostic methods or new tracers, we understand the importance of providing the very best care.

Intelligent technology

Our most advanced PET in a premium MR

At the heart of SIGNA™ PET/MR AIR™, you'll find the industry's most sensitive Time of Flight (ToF) PET detector and an exceptional MR-based attenuation correction (MRAC) with GE HealthCare's unique Zero-Echo-Time (ZTE) acquisition. This means you might see the smallest lesions and quantitate more accurately for each patient.



This advanced hardware runs on our latest MR 30 for SIGNA™ release, available across all of our SIGNA™ MR scanners. It features AIR™ Recon DL, our pioneering deep-learning-based reconstruction algorithm that enables pin-sharp image quality while simultaneously reducing scan time.

